

**PRODUCT: ERYTHROMYCIN ETHYLSUCCINATE 10% GRANULES/POWDER  
(Taste Masking Granules/Powder)**

**COMPOSITIONS:**

Each 100 gm Taste Masked Granules Contains:

Erythromycin Ethylsuccinate USP  
Eq. to Erythromycin ----- 10 gm  
Excipient Q.S

**SPECIFICATIONS:**

S. No.	TEST	SPECIFICATION
1	Description	White to slightly off white free flowing granules, which after reconstitution gives a white to off white coloured & flavoured suspension.
2	Identification by TLC	Should be comply
3	Loss on drying at 60 °C for 3 hours.	Limit: Not more than 5.0%
4	pH (Reconstituted suspension )	Limit : 7.0 to 9.0
5	Assay: Each 100 gm Contains Erythromycin Ethylsuccinate USP Eq. to Erythromycin .....10 gm	9.0 gm to 12.0 gm (90.0% to 120.0%)

**NOTE: Customized specification & flavour will be available on request.**